

Cleaning Your Spacer



1. Pull off protective cap and remove the end cap.
2. Rotate the mouthpiece to disassemble.
3. Once disassembled, gently sweep the spacer and mouthpiece through warm soapy water for 1 minute and allow both to soak for 10 minutes.
4. Repeat step 3 in clean warm water and shake out excess moisture.
5. Let all parts air dry thoroughly prior to reassembling.
6. To reassemble, place mouthpiece on chamber and rotate to lock.

Please contact Dr. Boas or a member of his clinical staff at 847-998-3434 for additional assistance.



CHICAGO
PULMONARY SPECIALISTS
THE SCIENCE OF BREATHING