

J. G. ROSENBAUM POCKET VISION SCREENER

95

distance
equivalent
20
/

874

2843

Point
Jaeger
26 16
20
/

638 E W E X O O

14 10
20
/

8745 E M W O X O

10 7
20
/

63925 M E E X O X

8 5
20
/

428365 W E M O X O

6 3
20
/

374258 W W X X O

5 2
20
/

8 X X X X W W X O O

4 1
20
/

4 J X T I I W W X O O

3 1+
20
/

Card is held in good light 14 inches from eye. Record vision for each eye separately with and without glasses. Presbyopic patients should read thru bifocal segment. Check myopes with glasses only.